


ADVANCED ENDOBRONCHIAL SOLUTIONS

VOLUTAM™ Endobronchial Coil System



NON-BLOCKING BRONCHOSCOPIC LUNG VOLUME REDUCTION

Clinical product overview for healthcare professionals

Appropriate Patient Selection

Data-driven. Patient-focused.

Volutam coil therapy may be considered for carefully selected patients with severe emphysema and marked hyperinflation. Patient evaluation should integrate symptoms, pulmonary physiology, CT/QCT findings, comorbidities and multidisciplinary risk-benefit assessment. Treatment suitability and technique must follow the current Volutam IFU.


CLINICAL ASSESSMENT

Optimized pharmacologic therapy
Pulmonary rehabilitation
Persistent symptoms / exercise limitation
FEV1, RV, TLC, RV/TLC
Gas exchange / vascular status
Cardiac and overall comorbidity

QCT-LED PLANNING

Lobar volume - anatomical burden
LAA-950 / LAA-910 - emphysema burden
Heterogeneity - target vs non-target
Fissure map - context, not a coil gate
Follow-up QCT - redistribution / safety

WHEN COIL CONSIDERATION BECOMES RELEVANT



WARNING - USE CAUTION / AVOID

Large focal or giant bullae; paraseptal-predominant emphysema; airway-predominant disease; clinically significant pulmonary hypertension; major uncontrolled cardiac disease; active infection or unstable exacerbation. Exact criteria and thresholds must follow the current Volutam IFU.

DECISION PRINCIPLE

No single CT or physiologic variable should be used in isolation. Selection is a multidisciplinary clinical decision integrating symptoms, physiology, imaging, comorbidity and the current approved labeling.

Hyperinflation



From air trapping to respiratory limitation

In emphysema, destruction of alveolar walls reduces elastic recoil. Expiratory flow limitation promotes air trapping, progressive hyperinflation and diaphragmatic flattening. The result is a rising mechanical burden on breathing and reduced exercise capacity.

HYPERINFLATED LUNG

TARGETED VOLUME REDUCTION

THE CYCLE OF OVERINFLATION




1 OVERDISTENSION
Alveoli lose elasticity and remain overdistended.



2 AIR TRAPPING
Expiratory flow limitation leaves air behind.



3 HYPERINFLATION
Lung volumes rise; the diaphragm flattens.



4 REDUCED FUNCTION
Mechanical efficiency and exercise capacity decline.

MECHANICAL CONSEQUENCE

Progressive air trapping increases lung volume, flattens the diaphragm and reduces mechanical efficiency. Targeted parenchymal gathering is intended to reduce the burden of hyperinflation in appropriately selected patients.

VOLUTAM™ CLINICAL OBJECTIVE

Mechanical gathering of targeted emphysematous parenchyma, with the aim of reducing hyperinflation and improving respiratory mechanics. Treatment suitability and technique must follow the current IFU.

Targeted volume reduction

Non-blocking. Mechanical. Image-guided.

Unlike one-way valves, coils do not aim to occlude an entire lobe. The deployed nitinol implant recovers its memorized form, gathers diseased parenchyma and increases local tissue tension. The mechanism therefore does not require complete lobar blockage.

1 PLAN

Physiology +
CT/QCT
target
selection

2 NAVIGATE

Bronchoscopy
+ fluoroscopy

3 DELIVER


Advance the
straightened
implant

4 DEPLOY

Release; coil
recovers
pre-set form

5 ASSESS

Confirm
position and
procedural
safety



CLINICAL DISTINCTION

Collateral ventilation is central to valve eligibility; coil suitability still requires comprehensive imaging, physiology and multidisciplinary assessment. Treatment suitability and technique must follow the current Volutam IFU.

Published coil-system evidence

Benefits and risks shown together

INDIVIDUAL PARTICIPANT DATA META-ANALYSIS • 680 PATIENTS

+0.07 L FEV 1 6 months	-0.36 L RV 12 months	-9.8 pts SGRQ 12 months	+38 m 6MWD 3 months
--	--	---	---

RENEW RCT • 315 PATIENTS • 12 MONTHS


+14.6 m 6MWD	-8.9 pts SGRQ	+7.0% FEV 1	34.8% vs 19.1% Major complications
--	---	---	--

INTERPRETATION

These values summarize published coil-system evidence and should not be presented as Volutam-specific efficacy data unless supported by an appropriate bridging rationale and regulatory review. Major complications were more frequent with coil treatment than usual care in RENEW.

The implant

Shape-memory nitinol for non-blocking bronchoscopic lung volume reduction.



SHAPE-MEMORY NITINOL

Delivered straightened; designed to recover its pre-set three-dimensional form after deployment.

NON-BLOCKING CONCEPT

Mechanical tissue gathering rather than complete lobar occlusion.

TARGETED PLACEMENT

Selected subsegmental airways under bronchoscopic and fluoroscopic guidance.

IMPLANT LENGTHS

100 mm

125 mm

150 mm

Product configurations and availability may vary by market. Refer to the current IFU and ordering information.

Precision delivery

Dedicated system for controlled implantation.


- 1

FORCEPS
Purpose-designed for coil grasping and positioning
- 2

GUIDE WIRE
Specially designed, nitinol wire-supported
- 3

APPLICATION CATHETER
For coil delivery
- 4

CARTRIDGE
Sterile packaging and protection



FORCEPS O.D. Ø 1.8 mm Opening width 5.5 mm Working length 1600 mm Working channel $\geq$ Ø 2.0 mm PURPOSE-DESIGNED Specially formed jaw geometry is designed to grasp the Volutam coil during implantation, supporting controlled positioning and adjustment.	GUIDE WIRE Overall length 1250 mm Outer diameter 1.7 mm Outer material 316LVM stainless steel Inner support Nitinol wire SPECIAL DESIGN 316LVM stainless-steel outer construction with internal nitinol wire support, engineered as a dedicated guide wire for the delivery system.	APPLICATION CATHETER Shaft length 100 cm Inner diameter 2.0 mm (0.081 in) Technical specifications as provided for the delivery components.
---	---	---

SINGLE-USE COMPONENTS
Supplied sterile | For use in bronchoscopic procedures

Refer to the current IFU for final specifications, nomenclature and instructions for use.

Risk management

Safety belongs in the main clinical history.

Clinical Risks and Considerations

- HEMOPTYSIS**

Clinically significant bleeding can occur.
- COPD EXACERBATION**

May occur early post-procedure.
- PNEUMONIA / LRTI**

Prominent signal in randomized trials.
- PNEUMOTHORAX**

Recognized post-procedural risk.
- COIL-ASSOCIATED OPACITY**

Radiographic phenomenon that may mimic infection.

MR CONDITIONAL



- 1.5 T / 3 T**

Field strength
- ≤ 2 W/kg**

Whole Body SAR
- Normal Mode**

Operating mode
- 15 min**

Maximum active scan

Refer to the current approved MRI labeling, implant card and IFU for complete MR Conditional conditions and safety information.

Clinical Evidence

Selected references supporting the clinical context presented in this catalog.

1. Sciurba FC, et al. RENEW. JAMA. 2016;315:2178-2189.
2. Deslée G, et al. REVOLENS. JAMA. 2016;315:175-184.
3. Shah PL, et al. RESET. Lancet Respir Med. 2013.
4. Hartman JE, et al. Systematic review and individual participant data meta-analysis. Respiration. 2022.
5. Slebos DJ, et al. Best Practice Recommendations. Respiration. 2018.
6. van Geffen WH, et al. Surgical and bronchoscopic pulmonary function-improving procedures in lung emphysema. Eur Respir Rev. 2023.
7. SEPAR Clinical Protocol on Endoscopic Lung Volume Reduction for Severe Emphysema. Open Respir Arch. 2026.



For healthcare professionals. Product availability, approved indications and labeling vary by market.
Refer to the current Instructions for Use.